

DIGITAL ASSESSMENT KIT (DAK)

Life-saving diagnostic equipment every yacht needs

Better assessment leads to better assistance.

MedAire's Digital Assessment Kit offers the most important diagnostic tools your crewmembers need to assess and mitigate medical situations, sending readings directly to MedLink.



KIT CONTENTS

BLOOD PRESSURE MONITOR Small and lightweight, easily adjustable to different sizes of people.

12 LEAD ECG RECORDER approved by the FDA & CE. It is easy-to-use, requiring minimal training. It is not a harness and can be easily applied to the chest of the guest/crewmember.

DIGITAL GLUCOMETER Less-invasive device which provides accurate results and requires no calibration

PULSE OXIMETER Provides instant results by clipping onto a person's fingertip to detect oxygen saturation, blood perfusion and respiratory rate.

CONTACTLESS/INFRARED THERMOMETER Takes temperature readings without physical patient contact.

IPAD With MedAire's New App, providing common yachting medical scenarios to assist the responder.

WHY CHOOSE THE DAK?

The most significant differentiator is the easy-to-use 12-Lead ECG recorder. Most devices of this size offer 1-Lead ECG which cannot supply enough information for a doctor to accurately assess the condition of the heart.

SPECIFICATIONS

LENGTH	203 mm	8 in
HEIGHT	102 mm	4 in
WIDTH	305 mm	12 in
WEIGHT	1.8 kg	4 lb


DIGITAL ASSESSMENT KIT (DAK)

HOW DOES THE DAK WORK?

The DAK is a tried and tested device, built by MedAire from the learnings of thousands of our assistance cases. It is already used by many of our BGA (private jet) clients.

Using MedAire's brand new MediConnect app, designed by our medical experts, especially for superyachts, responders can choose from a variety of scenarios from chest pain to vomiting and diving injuries, answering questions specific to the scenario.

Following the prompts from the App, responders can then transmit any vital signs readings to MedAire for diagnosis and evaluation and even when WiFi is not available, MediConnect can still be used to gather vital signs data and receive advice to improve patient care.



1. Activate the DAK and choose the scenario

2. Provide immediate care based on prompt

3. Provide any additional information requested

4. Transmit readings and call MedLink

DAK E-LEARNING

Our 1 hour e-learning is an interactive 6 module course. It is designed for crewmembers who have a basic knowledge of MedAire's DAK, helping them to review the components of the kit, learn how to look after it and how to operate it.

We recommend completing the training every year and in addition our instructors can incorporate DAK scenarios into your onboard training.

REPLENISHMENT



Expired items will be added to your next resupply:
- Pads: 36 months, if unused
- Test strips, new batteries & lancets: every 12 months.

REQUEST A QUOTE

For more information on pricing, reach out to your Account Manager.